

[HemoRO 3000]

Reliable Small Reverse Osmosis System



FUNCTION

The Small Reverse Osmosis Unit **HemoRO 3000** satisfies in home and acute dialysis for many years.

The **HemoRO 3000** reliable operating in many dialysis centres all over the world is even better and safer now.

Many suggestions of our customers, e.g. improved mobility, optimized noise level or higher permeate performance (50 or 70 l/h) were included in the development.

Your advantages at a glance:

- ▣ High microbiological safety through ring main to dialysis machine, semi-automatic chemical disinfection and fully automatic standby flushing
- ▣ Flexible use through hard water operation capability up to 25 °dH (4,46 mmol/l)
- ▣ Safe handling through colour coded tubing
- ▣ Easy handling through clear 2-button user interface
- ▣ Integrated pre-filter with pressure reducer

EQUIPMENT

The optional trolley enables an ergonomic transfer to the treatment place. Furthermore you can store the cables and tubing for better organisation and space-saving.

The DWA sampling set ensures an easy and standardised microbiological sampling.

For decalcification we advise the use of the DWA citric acid injection set.



DWA
Pollet Medical Group

DWA GmbH & Co. KG
The Dialysis Water Specialist

Großer Sand 8
76698 Ubstadt-Weiher

Phone: +49 7251 6900-0

info@dwa-online.com
www.dwa-online.com