

[Water Guard]

Safe – Cost Saving – Easy to Maintain



FUNCTION

The **DWA Water Guard** was developed especially for the use in the medical water processing. It offers optimal protection against unintended water leakage in your osmosis room – **even by power failures!**

If the water sensor comes into contact with liquids, the motor rotary valve closes within less than 10 seconds and interrupts the water supply for the reverse osmosis and all previous components like e.g. softeners.

With power failures, the ball valve is securely closed through an integrated “power pack” (not a storage battery, not a battery). As opposed to standard processes, the function of this DWA technology cannot be blocked by lime or other deposits.

As soon as the mains voltage is available again, the motor rotary valve automatically opens if no water leakage has been recognised.

With a voltage or control failure, the motor rotary valve can also be opened manually.

Your advantages at a glance:

- Optimal building protection even by power failures

If you have questions about the **DWA Water Guard** please feel free to contact our Export Department.



DWA
Pollet Medical Group

DWA GmbH & Co. KG
The Dialysis Water Specialist

Großer Sand 8
76698 Ubstadt-Weiher

Phone: +49 7251 6900-0

info@dwa-online.com
www.dwa-online.com